

Universal Multi Rack



Product Categories: [Clamps, Holders Racks Stands](#), [Laboratory Plastic Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/universal-multirack/>

Product Description

Universal Multi Rack - Laboratory Plastic ware Equipment

This rack is an excellent lab apparatus compatible with different sizes of tubes. Four sides of this rack can be used for tubes of 30 mm, 20 mm, 17 mm & 12 mm.

Material : Polypropylene

High grade Polypropylene is used for its manufacturing.