

# Universal Oven



**Product Categories:** [Hot-Air Ovens](#), [Lab Heating- Equipment](#), [Laboratory Equipments](#), [Ovens](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/universal-oven/>

## Product Summary

Optional Accessories

- Timer - Air Circulating Fan - Digital Temperature Indicator - Digital Controlled- Cum- Indicator

## Product Description

Universal Oven - Scientific Equipment

Temperature up to 250 C. Heating elements are placed in ribs at the bottom and sides. Double walled inside anodized Aluminium/stainless steels. the door has synthetic rubber gasket. outside Mild steel Painted in epoxy powder coating with perforated adjusted shelves two indicators and on/off switch. To work on 220/230 volts A.C with on number Air Circulation fan Digital Temperature Controller Cum Indicator.

Riding on unfathomable volumes of industrial proficiency, we are presenting a wide compilation of UNIVERSAL OVEN Provided optimum quality products are made accessible in several sizes.

Presented products are precisely made-up using the supreme quality material and contemporary technology under the command of competent personnel.

### Features

Precise design

Superior finish

Light weight