

Universal Vibration Apparatus



Product Categories: [Engineering Equipment](#), [Theory Of Machine Lab](#)

Product Page: <https://www.labappara.com/product/universal-vibration-apparatus/>

Product Description

Universal Vibration Apparatus

The apparatus is provided with comprehensive unit to perform the vibration experiments. A universal frame is provided upon which quick and easy assembly of various experiments can be done. The unit is self-contained to safely store spares. The students can easily assemble the experiments and study the theory of vibrations practically.

Technical Details:

Exciter Unit: With FHP DC motor with speed control facility

RPM Measurement: Digital RPM Indicator with proximity sensor

Ordinary Chart Recorder: For recording frequency and amplitude of vibration

Stop Watch: Electronic Stop Watch

Instruction Manual: An English instruction manual will be provided along with the apparatus