

Vacuum Pump



Product Categories: [Pathology Lab](#)

Product Tags: [Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/vacuum-pump/>

Product Description

Vacuum Pump - Laboratory Equipment

A pump used for creating a vacuum.

Vacuum pump is consists of single stage, oil sealed, air cooled.

V- Belt driven and mounted on base plate with belt guard.

Vacuum pump with motor 220volt AC single phase.

Vacuum displacement 32 litre / Min.

Vacuum pump is consisting of double stage 75 litre cap.

Applications:

Chemical Industry

Vacuum Testing

Vacuum Lifting

Distillation

Filtration

Vacuum Costing

Salient Features:

Compact Design

Low Noise Level

Maintenance Free

Easy to handle

Continuous Duty

Totally Dry Compressed Air

Low Running Cost