

Vapour Pressure Capillary Type Dial Thermometer



Product Categories: [Chemistry Equipments](#), [Thermometer & Hygrometer](#)

Product Tags: [chemistry lab instruments](#), [physics lab instruments](#)

Product Page:

<https://www.labappara.com/product/vapour-pressure-capillary-type-dial-thermometer/>

Product Description

Vapour Pressure Capillary Type Dial Thermometer - Laboratory Equipment

100 mm dials size x 2

Meter PVC Capillary,

Range- 0 to 500C *

Range -0 to 6000C *

Features:

Durability

Excellent performance

Highly accuracy

Technical Specifications :

Dial size

100, 150 mm

Case

W.P Cast AL. or Stainless Steel

Mounting

Surface, Panel, bottom or back entry

Type

Capillary distance or Rigid Stem

Capillary

3 Mtr. PVC, Copper or stainless Steel

Bulb, Stem & Connection

Material

Brass or Stainless Steel

Bulb & Stem Length

Std. 12" long or as required

Bulb & Stem Dia

1/2" OD

Connection

1/2" BSP (M) or as required