

Venturi Flume



Product Categories: [Engineering Equipment](#), [Hydraulic Lab](#)
Product Page: <https://www.labappara.com/product/venturi-flume/>

Product Description

Venturi Flume

Specification :-

Venturi flume for the experimental flume

Venturi flume consisting of 1 base plate, 2 side elements, 1 clamping device

Side elements with sealing lips

Technical Data Venturi flume

L x W x H : 1000 x 304 x 450 mm

Narrowest cross-section, W x H: 152 x 430 mm

Side element

L x W x H : 828 x 76 x 430 mm