

Vertical Pipette Stand



Product Categories: [Laboratory Plastic Ware](#), [Pippetes Stand](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/vertical-pipette-stand/>

Product Description

Vertical Pipette Stand - Scientific Laboratory Plastic ware

Pipette Stand (Vertical):- This is a stable, non-corrosive, detachable stand which can hold 28 pipettes with maximum dia 14 mm. The heavy base measures 225 mm in diameter and a metal reinforced square Polypropylene rod measuring 375 mm in length is screwed in the base. The top disc that can be fixed on the rod has holes for holding pipettes while the lower one has pilot holes where the pipette tips rest and liquids drain down to the base. These disc can be adjusted on the rod according to the length of the pipette. The arrows provided on the disc should be kept parallel to ensure that corresponding holes are being centered.

Particulars	Packing
For 28 Pipettes	1 Pcs