

## Vicat Apparatus



**Product Categories:** [Pharmaceutical Equipments](#)

**Product Tags:** [Laboratory Equipment](#)

**Product Page:** <https://www.labappara.com/product/vi-cat-apparatus/>

### Product Description

#### Vicat Apparatus - Laboratory Equipment

A Vicat apparatus was used to measure standard consistence, and initial and final setting time after the reactants were mixed together and placed in a mould. The apparatus is provided with dash-pot arrangement for jerk free lowering of needles

#### **Available with**

Vicat mould 60 mm top dia, 70 mm bottom dia and 40 mm height

Vicat Mould 70 mm top dia, 80 mm bottom dia and 40 mm height

#### **Features:**

Rigid structure

Impeccable finish

Durability