

Vicat Apparatus



Product Categories: [Pharmaceutical Equipments](#)

Product Tags: [Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/vi-cat-apparatus/>

Product Description

Vicat Apparatus - Laboratory Equipment

A Vicat apparatus was used to measure standard consistence, and initial and final setting time after the reactants were mixed together and placed in a mould. The apparatus is provided with dash-pot arrangement for jerk free lowering of needles

Available with

Vicat mould 60 mm top dia, 70 mm bottom dia and 40 mm height

Vicat Mould 70 mm top dia, 80 mm bottom dia and 40 mm height

Features:

Rigid structure

Impeccable finish

Durability