


Vibration Magnetometer (wood box)



Product Categories: [Educational equipments](#), [Electromagnetism](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/vibration-magnetometer-wood-box/>

Product Description

Vibration Magnetometer (wood box) - Physics Instrument

Consisting of a superior quality polish wooden box slotted, and carries a glass tube fitted with simple torsion head to support, by means of a silk fiber, a strip holding magnet etc. Complete with one bar magnet and one brass bar. All metal parts are of magnet and one brass bar. All metal parts are manufactured from superior quality brass alloy.

With Aluminum Fittings - Consisting of a superior quality polish wooden box slotted, and carries a glass tube fitted with simple torsion head to support, by means of a silk fiber, a strip holding magnet etc. Complete with one bar magnet and one brass bar. All metal parts are of magnet and one aluminum bar. All metal parts are manufactured from superior quality Aluminum alloy.

With Brass Fittings - - Consisting of a superior quality polish wooden box slotted, and carries a glass tube fitted with simple torsion head to support, by means of a silk fiber, a strip holding magnet etc. Complete with one bar magnet and one brass bar. All metal parts are of magnet and one brass bar. All metal parts are manufactured from superior quality brass alloy.