

## Voltage Regulation Using IC 317



**Product Categories:** [Characteristics And Application](#), [Characteristics And Application](#), [Electronics](#), [Engineering Equipment](#)

**Product Page:**

<https://www.labappara.com/product/voltage-regulation-using-ic-317/>

### Product Description

#### Voltage Regulation Using IC 317

##### **Features :**

Instrument comprises of AC Power Supply, three meters to measure Output Voltage, Output Current & AC Ripples directly, two PN Junction Diodes, Filter, Regulation circuit using IC 317, load Resistances selectable using band switch & all important connections brought out on 4mm Sockets.