

Volume Relationship Set



Product Categories: [Geometry Items](#), [Math Lab](#)

Product Tags: [mathematics lab equipments](#)

Product Page: <https://www.labappara.com/product/volume-relationship-set/>

Product Description

Volume Relationship Set - Mathematics Lab Equipments

To teach volume relationship between 3D solids, this Set consists with 10 cm dimension 3D solids that allow to fill or any dry material to demonstrate volume relationship. This set consists of cone cylinder, Square prism and pyramid and a sphere equal to the inner diameter of cylinder or square prism.