

Volumetric Flask Borosilicate Glass



Product Categories: [Laboratory Glassware](#)

Product Tags: [glassware](#)

Product Page: <https://www.labappara.com/product/volumetric-flask-borosilicate-glass/>

Product Description

Volumetric Flask Borosilicate Glass - Laboratory Glassware Equipment

Made from **Borosilicate Glass**, with heavy tubing neck. Calibrated to contain, with fitted polyethylene stopper.

Feature -

- Capacity is 500ML - Borosilicate Glass Flask - Polyethylene Stopper