

## Volumetric Flask



**Product Categories:** [Laboratory Plastic Ware](#), [Liquid-Measurement](#)

**Product Tags:** [Laboratory Equipment](#), [Laboratory plastic ware](#)

**Product Page:** <https://www.labappara.com/product/volumetric-flask/>

### Product Description

Volumetric Flask - Scientific Laboratory Equipment

Volumetric Flasks, blown in Polypropylene, are strong, unbreakable and have good contact clarity. These autoclavable flasks are provided with screw caps to prevent leakage.

Volumetric Flasks are blown in Polypropylene and are strong, unbreakable and have good contact clarity.

These autoclave Volumetric Flasks are provided with screw caps to prevent leakage.

These Volumetric Flasks are economic as well.

#### **Features:**

Easy to operate

Lightweight

Accurate measurement

**Material:** Polypropylene

**Available in:** 25 ml, 50 ml, 100 ml, 250 ml, 500 ml, 1000 ml