

W.B.C / R.B.C Pipettes



Product Categories: [Pathology Lab](#)

Product Tags: [Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/w-b-c-r-b-c-pipettes/>

Product Description

W.B.C / R.B.C Pipettes - Laboratory Equipment

A pipettes glass tube having a uniform bulb which contains a red bead (for R.B.C) and a white bead (for W.B.C).

The tube is white back with permanent graduations and the tip of the tube is finely ground and beveled to ensure precise sampling.

The other end of the glass tube is fitted with a 20 cm latex or transparent PVC mouth piece tubing.

Features:

High quality

Fine finish

Longer service life

High performance

Precise construction