

Water Bath copper



Product Categories: [Laboratory Equipments](#), [Water Bath& Stills](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/water-bath-copper/>

Product Description

Water Bath copper - Scientific Laboratory Equipment

Hemispherical, made of thick copper sheet, with a set of concentric rings and 2 handles.

Hemispherical made from thick sheet with side rings bath with side strips and concentric rings for heating beaker.

Water baths can hold often temperatures within a tenth of a degree Celsius.

Material: Copper Weight: 3 to 4 kg

Size-4" Size-5" Size-6" Size-8