

Water Bath (Electric)



Product Categories: [Laboratory Equipments](#), [Water Bath& Stills](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [Manufacturer of Laboratory equipment](#), [math lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/water-bath-electric/>

Product Description

Bath (Tissue Floation) - Scientific Laboratory Equipment

This water bath is thermostatically controlled. Constructed with double stainless steel walls and inner chamber made of black oxidized aluminum sheet, this device is perfect for measuring tissue culture and enzyme reactions. Temperature range of this device is from ambient to 90 degrees C. It works on 220 V AC supply and the front panel has temperature regulation knob along with power switch and pilot lamp.

200 x 75 mm

250 x 125 mm

Features:

Compact size

Robust construction

Excellent functioning

Easy operations

Superior functioning

Superior performance