

Water Bath Rectangular (Electric single 6 hole)



Product Categories: [Laboratory Equipments](#), [Water Bath& Stills](#)

Product Page:

<https://www.labappara.com/product/water-bath-rectangular-electric-single-6-hole/>

Product Description

Water Bath Rectangular (Electric single 6 hole) - Scientific

Laboratory Equipment

Thermostatically controlled. Double walled. Glass window on tow opposite sides thermometer with relay, complete with stirrer. High precision water bath.

Temperature range from ambient temperature to 100 C. Temperature controlled by Mercury contact thermometer with a sensitivity of + 0.1 0 C. Suitable for all laboratory applications requiring constant temperature with high degree accuracy. Double walled inside Stainless steel and outside mild steel sheet. A glass window is fixed on the front. Complete with stirrer, with speed regulator.

Made of Stainless steel sheet. Top of the bath and concentric rings are also of Stainless Steel. Heating elements have an automatic cut off device when the bath gets dry. Electrically operated on 220/230 volts A.C. with arrangement of 6 holes