

Water Taps



Product Categories: [Chemistry Equipments](#), [Lab Supplies](#)

Product Tags: [chemistry lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/water-taps/>

Product Description

Water Taps - Scientific Equipment

We are engaged in offering Water Tap , which are either made of brass or stainless steel and are used in hospitals, colleges, laboratories and research institute. Available in black oxidized. Made of brass, inlet tail is threaded ½ inch Bsp outlet tube is made of 15 mm, brass tubing which is coupled to the tap standard by means of a sleeve not allowing the neck to swivel or be locked in a desired position. Leak proof operation, bench mounting.

Laboratory water taps.

Excellent quality.

Precision engineering.

Single Swan Goose Neck Water Valve are made from forging grade material to avoid Porosity and Leakage.

Feature-

Swan Neck type.

Gun Metal C.P leak proof.

Provide with One way, Two way, & Three way as demanded by customer.