

Water Voltmeter



Product Categories: [Physiology Apparatus](#)

Product Tags: [Scientific Equipment](#)

Product Page: <https://www.labappara.com/product/water-voltmeter/>

Product Description

Water Voltmeter - Scientific Equipment

Cup shaped with fittings and two graduated test tubes. The glass vessel is mounted on a sturdy

wooden base fitted with terminals and eureka wire electrodes. Supplied with brass metal test tube

holder and graduated tubes.

A voltmeter is a scientific instrument used for measuring quantity of electricity (electric charge) through electrolytic action. The SI unit of quantity of electricity is the coulomb. The voltmeter should not be confused with a voltmeter, which measures electric potential.

Water voltmeter in with glass tube and with bakelite part also.

Water voltmeter -do- with platinum electrode.

Features:

Easy to use

Long working life

Precisely engineered