

Wave Motion Apparatus



Product Categories: [Physics Equipments](#), [Sound And Waves Oscillations](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/wave-motion-apparatus/>

Product Description

Wave Motion Apparatus- Physics Lab Equipment

Wave Motion Apparatus Is used In Schools and colleges for visualizing the Propagation of wave.

Wave Motion Apparatus is used for demonstrating the movement of progressive transverse waves.

The Wave Motion Apparatus comprises of small balls fitted on a wire rods.

Wave Motion Apparatus which is ideally used for demonstrating longitudinal & transverse motion.

It consists of 18 eccentric discs supporting a series of metal rods.

On revolving the handle, transverse waves are obtained with bent.

They are straight rods and complete on wooden base.

Feature

Wave Motion Apparatus with plastic discs fitted.

Nicely polished wooden case

Provide with 18 & 24 pulleys.

Robustness

Perfect finish

Precisely designed