

Wide Mouth square Bottle



Product Categories: [Bottles And Asssesories](#), [Laboratory Plastic Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/wide-mouth-square-bottle/>

Product Description

Wide Mouth square Bottle - Scientific Laboratory Equipment

Material : Polypropylene

This autoclavable bottle, made of polypropylene, can be used for storing & shipping chemicals & reagents. Being wide mouthed, it is always convenient to pour in & take out the reagents. The built in depressions on two sides provide excellent grip while holding this bottle.

A bottle is a rigid container with a neck that is narrower than the body and a mouth. By contrast, a jar has a relatively large mouth or opening which may be as wide as the overall container.

Particulars	Packing
4000 ml	1 Pc