

Wide Mouth Wash Bottle



Product Categories: [Bottles And Assesories](#), [Laboratory Plastic Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/wide-mouth-wash-bottle/>

Product Description

Wide Mouth Wash Bottle - Scientific Laboratory Equipment

A bottle is a rigid container with a neck that is narrower than the body, and a "mouth".

Poly-lab Wash Bottle is made of Low Density Polyethylene, which imparts a Translucent and unbreakable character. Solutions can be poured easily by squeezing the Wash Bottle lightly. The Cap is fitted with a Flexible Polyethylene Delivery Tube which can be aimed wherever required.

Material : Polyethylene

Particulars	Packing
125 ml	12 Pcs
250 ml	12 Pcs
500 ml	6 Pcs
1000 ml	6 Pcs