

Wintrobe Tube



Product Categories: [Pathology Lab](#)

Product Tags: [Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/wintrobe-tube/>

Product Description

Wintrobe Tube - Scientific Equipment

It is a common Wintrobe tube, to perform the test, anti-coagulated blood was traditionally placed in an upright tube, known as a Wintrobe tube, and the rate at which the red blood cells fall was measured and reported in mm/h.

These tubes can be centrifuged.

The quantity of fluoride-oxalate solution per tube is very critical.

Little less quantity will result in clotting of blood sample and little high quantity will result in clumping of blood cells.

Features:

Easy to use

Eco-friendly

Recyclable

Perfectly designed

Durable