

## Wintrobe Tube



**Product Categories:** [Pathology Lab](#)

**Product Tags:** [Laboratory Equipment](#)

**Product Page:** <https://www.labappara.com/product/wintrobe-tube/>

### Product Description

#### Wintrobe Tube - Scientific Equipment

It is a common Wintrobe tube, to perform the test, anti-coagulated blood was traditionally placed in an upright tube, known as a Wintrobe tube, and the rate at which the red blood cells fall was measured and reported in mm/h.

These tubes can be centrifuged.

The quantity of fluoride-oxalate solution per tube is very critical.

Little less quantity will result in clotting of blood sample and little high quantity will result in clumping of blood cells.

#### **Features:**

Easy to use

Eco-friendly

Recyclable

Perfectly designed

Durable