


Wire Pattern Test Tube Stand



Product Categories: [Clamps, Holders Racks Stands](#), [Laboratory Plastic Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/wire-pattern-test-tube-stand/>

Product Description

Wire Pattern Test Tube Stand - Scientific Plastic ware Equipment

Test tube stand is a laboratory equipment that is used to hold upright multiple test tubes at the same time. It is most commonly used when various different solutions are needed to work with simultaneously, for safety reasons, for safe storage of test tubes, and to ease the transport of multiple tubes.

For cleaning purposes Test Tube Stand can easily be disassembled. The base is in the form of a grid with square openings and has fold up sides and bars that lock into rack tops.

Particulars	Packing
13 mm X 36 Tubes	4 Pcs
16 mm X 36 Tubes	4 Pcs