


Sodium Vapour Lamp Wooden Box



Product Categories: [Educational equipments](#), [Electronics](#)

Product Tags: [electronics](#), [optics](#)

Product Page:

<https://www.labappara.com/product/wooden-box-sodium-vapour-lamp/>

Product Description

Sodium Vapour Lamp Wooden Box - Laboratory Equipment

The box is rectangular in shape and provided with slits on three sides equally spaced to enable three experiments to be performed simultaneously using one light source

Each slit is provided with a cover so that it can be closed when not required

Bulb holders of finest quality are fitted in wooden box with special cord Quality:

Reliable

Best quality

Effective service

Available with 35 watt bulb and transformer Transformer:

High grade laminations are used to avoid losses

The box is made from thick metallic sheet

An anchor make piano type switch is provided for on / off and a 1.6 m, 3 care flying lead is provided for main connection and output is via a special Anchor make 3-pin

socket (3-pin socket is provided to avoid shock due to the leakage of the current)
This transformer is an auto leak type, giving a starting voltage 400 V which automatically reduces when the lamp has struck.

Available in Watts :

(a) 35 watts

(b)55 watts