

Work Function Of Diode (Richardson Law)



Product Categories: [E/M Apparatus & Valve Experiments](#), [E/M Apparatus & Valve Experiments](#), [Electronics](#), [Engineering Equipment](#)

Product Page:

<https://www.labappara.com/product/work-function-of-diode-richardson-law/>

Product Description

Work Function Of Diode (Richardson Law)

Features :

Instrument comprises of DC Regulated power supply of 0-6 V & 0-300 V. Two volt meters & two current meters are mounted on the front panel. Valve diagram is printed on the panel & connection are brought out on sockets/terminals.