

Y(Tube) Connectors



Product Categories: [Connectors And Stirrers](#), [Laboratory Plastic Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/ytube-connectors/>

Product Description

Y(Tube) Connectors - Scientific Laboratory Equipment

These connectors add to the existing range of connectors. The crest & valley structure provided at all the ends of the connectors give excellent leak proof grip to the tube.

Features:

Robust construction,

Durability,

Dimensional accuracy

Consistent performance.

Particulars	Packing
Y Connector-6 mm	36 Pcs
Y Connector-8 mm	36 Pcs
Y Connector-10 mm	36 Pcs